

Catalog	Name	Specificity	Host Cell	Growth/Thaw Medium	Vector	Size
B23602614	Cyno ILT2 Cell Line (CHO-K1)	Cyno ILT2	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602615	Cyno ILT4 Cell Line (CHO-K1)	Cyno ILT4	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602612	Cyno TNFR2 Cell Line (CHO-K1)	Cyno TNFR2	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602601	Human B7H3 Cell line (CHO-K1)	Human B7H3	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602616	Human 2Ig-B7H3 Cell Line (CHO-K1)	Human B7H3	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602622	Human 4Ig-B7H3-HA Cell Line (CHO-K1)	Human B7H3	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23144302	Human BCMA Cell line (HEK293)	Human BCMA	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B23602626	Human CCR10 Cell Line (HEK293)	Human CCR10	HEK293	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6
B21927402	Human CD137 Cell line (CHO-K1)	Human CD137	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22990201	Human CD19 Cell line (3T3)	Human CD19	3T3	90%DMEM+10%FBS+2.5ug/ml puro	Lentiviral	2*6E6
B22990401	Human CD19 Cell line (CHO-K1)	Human CD19	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22990301	Human CD19 Cell line (HEK293)	Human CD19	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B23147201	Human CD20 Cell line (HEK293)	Human CD20	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B22997201	Human CD20 Cell line (CHO-K1)	Human CD20	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B23011501	Human CD24 Cell line (CHO-K1)	Human CD24	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B21927403	Human CD28 Cell line (CHO-K1)	Human CD28	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22359201	Human CD33 Cell line (CHO-K1)	Human CD33	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23012101	Human CD38 Cell line (CHO-K1)	Human CD38	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22339902	Human CD70 Cell line (CHO-K1)	Human CD70	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22428701	Human CD89 Cell line (CHO-K1)	Human CD89	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22776901	Human CD8A & CD8B Cell line (CHO-K1)	Human CD8A & CD8B	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22468801	Human CD8B Cell line (CHO-K1)	Human CD8B	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602606	Human CLEC5A-Flag Cell Line (CHO-K1)	Human CLEC5A	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B22468802	Human CMET Cell line (CHO-K1)	Human CMET	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22339901	Human CTLA4 Cell line (CHO-K1)	Human CTLA4	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602603	Human DLL3(AA312-618) Cell line(CHO-K1)	Human DLL3	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602602	Human DLL3 Cell line (CHO-K1)	Human DLL3	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B23602617	Human EPHA2(AA382-922) Cell Line (CHO-K1)	Human FCRH5	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6

B23149201	Human FCRH5 Cell line (CHO-K1)	Human FCRH5	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6
B23602604	Human GPC3 Cell line (3T3)	Human GPC3	3T3	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6
B22381101	Human GPRC5D Cell line (HEK293)	Human GPRC5D	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B21630701	Human GPRC5D Cell line (CHO-K1)	Human GPRC5D	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B21583004	Human GUCY2C Cell line (CHO-K1)	Human GUCY2C	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602607	Human LINGO1-Flag Cell Line (CHO-K1)	Human LINGO1	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B22131101	Human LRRC15 Cell line (CHO-K1)	Human LRRC15	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22416901	Human MICA Cell line (CHO-K1)	Human MICA	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22468803	Human MSLN Cell line (CHO-K1)	Human MSLN	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602608	Human MUC13-Flag Cell Line (CHO-K1)	Human MUC13	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602623	Human Flag-MUC16-HA Cell Line (CHO-K1)	Human MUC16	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B21583003	Human NaPi2B Cell line (CHO-K1)	Human NaPi2B	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22375101	Human Napi2B Cell line (3T3)	Human Napi2B	3T3	90%DMEM+10%FBS	Lentiviral	2*6E6
B23142201	Human Nectin4 Cell line (CHO-K1)	Human Nectin4	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22719701	Human PD-1 Cell line (CHO-K1)	Human PD-1	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22301201	Human PD-L1 Cell line (CHO-K1)	Human PD-L1	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B22190801	Human ROR1 Cell line (CHO-K1)	Human ROR1	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602627	Human SLC4A7(AA601-1131) Cell Line (CHO-K1)	Human SLC4A7	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B23602628	Human SLC9A1 Cell Line (CHO-K1)	Human SLC9A1	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B22428601	Human TCR/CD3 Cell line (CHO-K1)	Human TCR/CD3	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B23147301	Human TIGIT Cell line (CHO-K1)	Human TIGIT	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22529801	Human TMPRSS4 Cell line (3T3)	Human TMPRSS4	3T3	90%DMEM+10%FBS+2.5ug/ml puro	Lentiviral	2*6E6
B22529802	Human TMPRSS4 Cell line (HEK293)	Human TMPRSS4	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B22529803	Human TMPRSS4 Cell line (CHO-K1)	Human TMPRSS4	CHO-K1	90%DMEM/F-12+10%FBS	Lentiviral	2*6E6
B23602611	Human TNFR2 Cell Line (CHO-K1)	Human TNFR2	CHOK1	DMEM/F12+10%FBS+1% Dual Antibiotics+10ug/ml puro	Lentiviral	2*6E6
B22845401	Human TPBG Cell line (CHO-K1)	Human TPBG	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22845402	Human TPBG Cell line (HEK293)	Human TPBG	HEK293	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B23602613	Human TROP2(AA88-323)-HA Cell Line (HEK293)	Human TROP2	HEK293	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6
B23602625	Human TROP2 R87A/T88A-HA Cell Line (CHO-K1)	Human TROP2	HEK293	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6

B23149301	Human Trop2 Cell line (HEK293)	Human Trop2	CHO-K1	90%DMEM+10%FBS+2ug/ml puro	Lentiviral	2*6E6
B23149401	Human Trop2 R87A/T88A Cell line (HEK293)	Human Trop2	HEK293	90%DMEM+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B22866101	Mouse CD20 Cell line (CHO-K1)	Mouse CD20	CHO-K1	90%DMEM-F12+10%FBS+1%Dual Antibiotics	Lentiviral	2*6E6
B23602629	Mouse SLC4A4(AA455-932)Cell Line (CHO-K1)	Mouse SLC4A4	HEK293	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6
B23602630	Mouse SLC9A1 Cell Line (CHO-K1)	Mouse SLC9A1	HEK293	DMEM+10%FBS+1% Dual Antibiotics+2ug/ml puro	Lentiviral	2*6E6